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LAVANID[®] - Wundgel / LAVANID[®] - Wundgel V+

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en INSTRUCTIONS FOR USE

LAVANID[®] - Wound Gel / LAVANID[®] - Wound Gel V+ Read the instructions for use completely before using the product. Use the product only for its intended purpose. Please note the meanings of the symbols at the end of the instructions for use.

Note for patients: This product is used by a doctor or other healthcare professional for the treatment of wounds. Patients should be educated by a doctor or other healthcare professional before use. The type and duration of use is determined by the condition of the wound.

If you notice any adverse effects, contact your doctor or healthcare professional.

DESCRIPTION AND INTENDED PURPOSE

LAVANID[®] Wound Gel and LAVANID[®] Wound Gel V+ are sterile, ready-to-use wound gels for cleaning and moistening wounds. LAVANID[®] Wound Gel V+ is more viscous than LAVANID[®] Wound Gel.

This product is suitable for acute wounds and for repeat-treatment of chronic wounds. Chronic wounds are often covered with cell debris and dead (necrotic) tissue and may be colonised by germs. Wound coatings can lead to delayed wound healing and should be removed during wound treatment, for example by debridement. Gentle wound cleansing and wound moistening is necessary as an accompanying measure. In individual cases, acute reactions of the immune system (anaphylactic reactions) have been reported. Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

ADVERSE EFFECTS

The ingredient polyhexanide can cause allergic reactions. In individual cases, acute reactions of the immune system (anaphylactic reactions) have been reported. Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE AND SHELF LIFE

The product must not be used after the expiry date. The expiry date is indicated on the label.

Temperature limit: Store between 2 °C and 30 °C.

Store in a dry place.

Do not use if the product is stored outside the recommended storage conditions.

Fragile, handle with care.

Keep out of reach of children.

Shelf life after opening: 1 month.

Due to the preservation, a self-life of one month is possible after the first opening of the container. A prerequisite for this is that the gel is protected from contamination. The container must be closed immediately after each use. Gel residue around the opening should be wiped off with a sterile compress. The container opening must be protected from contamination. The container must not come into contact with the wound. Containers which have come into contact with the wound or have been contaminated in any other way must be disposed of.

INFORMATION ABOUT OPENING THE CONTAINER

Information about opening the container can be found at the end of these instructions for use.

DISPOSAL INFORMATION

The disposal of unused products or waste material must be carried out in accordance with specific national regulations.

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PRECAUTIONS AND WARNINGS

STERNLE The product is sterile as long as the container is unopened and undamaged.

Do not use if the product shows signs of damage or if the packaging has been damaged or accidentally opened before use.

Do not sterilize/reprocess. This can lead to critical changes in performance, safety and handling of the product.

Do not use for injection or infusion.

Do not drink or swallow the product.

Do not mix with additives.

INTERACTIONS WITH OTHER PRODUCTS

When used together with other wound treatment products, the instructions for use of these products should be checked for possible interactions with the ingredients of the wound gel.

Due to incompatibilities with the ingredient polyhexanide, the wound gel must not be used together with anionic substances. These are, for example, substances that reduce the surface tension (anionic surfactants) or anionic gel formers or products containing PVP-iodine. In instances where there is prior contact, these substances must be rinsed liberally off of the wound.

ADVERSE EFFECTS

The ingredient polyhexanide can cause allergic reactions. In individual cases, acute reactions of the immune system (anaphylactic reactions) have been reported. Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE AND SHELF LIFE

The product must not be used after the expiry date. The expiry date is indicated on the label.

Temperature limit: Store between 2 °C and 30 °C.

Store in a dry place.

Do not use if the product is stored outside the recommended storage conditions.

Fragile, handle with care.

Keep out of reach of children.

Shelf life after opening: 1 month.

Due to the preservation, a self-life of one month is possible after the first opening of the container. A prerequisite for this is that the gel is protected from contamination. The container must be closed immediately after each use. Gel residue around the opening should be wiped off with a sterile compress. The container opening must be protected from contamination. The container must not come into contact with the wound. Containers which have come into contact with the wound or have been contaminated in any other way must be disposed of.

INFORMATION ABOUT OPENING THE CONTAINER

Information about opening the container can be found at the end of these instructions for use.

DISPOSAL INFORMATION

The disposal of unused products or waste material must be carried out in accordance with specific national regulations.

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