**NYLON SERALON® SUPRAMID**

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- **NYLON**
- **SERALON®**
- **SUPRAMID**

**DESCRIPTION**

NYLON, SERALON® and SUPRAMID are synthetic, nonabsorbable, sterile surgical suture materials of polyamide. Polyamides are polymers that are manufactured from amines and dicarboxylic acids or from amino-carboxylic acids or their functional derivatives. SUPRAMID in the cassette pack is supplied in an isopropanol-containing preserving solution.

Our polyamide suture materials are supplied in combinations with stainless steel needles or without needles. They comply with the underlying requirements of the Medical Device Directive 93/42/EEC and apart from marginally deviations with respect to the diameter (SUPRAMID), the current edition of the European Pharmacopoeia monograph (EP) for "Sutures, sterile non-absorbable (Fila non resorbilia sterilia)" or of the monograph "Nonabsorbable surgical suture" of the United States Pharmacopoeia (USP). In deviation from this standard, the maximum thread length in the case of the cassette packs is not 4, but 100 m.

- **NYLON** black

NYLON is a monofilament thread and is dyed with Hämatein (Logwood) black, Colour Index No. 75290. NYLON is available in very fine gauge sizes.

- **SERALON®** blue and undyed

SERALON® is a monofilament thread. As a result of its manufacturing process, it is especially pliable. SERALON® blue is dyed with Phthalocyaninblau, Colour Index No. 74160.

- **SUPRAMID** black and undyed

SUPRAMID is made of a twisted filament-bundle and a sheath (gauge sizes EP1 and thinner: monofilament) and therefore appears as a monofilament. SUPRAMID black is dyed with pigment black, Colour Index No. 77266.

**INDICATIONS**

Polyamide threads are suitable in principle for all wound sutures and ligatures in which nonabsorbable thread material is indicated and long-term stability

of the suture material is not required. Thanks to their pliability, low tendency to swell, and absence of a wick effect, SERALON® and SUPRAMID are among the most suitable of all suture materials for all types of skin suture. NYLON can be manufactured up to a gauge size of EP 0.1 and is therefore particularly suitable for sutures in microsurgery and ophthalmic surgery.

Polyamide threads can also be used as holding sutures and for marking purposes.

**USE**

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security (reg. SERALON®: see WARNINGS).

In order to prevent needle damage during suturing it is recommended that the needle be grasped at a point between the end of its first third (as seen from its end) and its midpoint.

**ACTION**

The function of a medical suture in wound care is to attach tissue to tissue or tissue to allogeneic material. Ligatures are a special type of suture intended to close hollow organs.

Polyamide threads are nonabsorbable. Nevertheless, threads left in tissues for a long time can lose tensile strength and eventually can fragment.

**CONTRAINDICATIONS**

Polyamide suture materials are not intended for use in the central vascular or nervous systems. Because of the loss of tensile strength that occurs in the body, polyamide threads must not be used when prolonged stability of the suture material is required.

**WARNINGS / PRECAUTIONS / INTERACTIONS**

The requirements of the suture material and the techniques required vary with the indication. In the choice of suture material, account should be taken of the in-vivo characteristics of the suture material concerned.

In view of the high elasticity of thicker monofilament polyamide threads, it is recommended that one to two additional knots be applied in order to achieve high knot-holding security; in the case of multifilament polyamide threads, a distance of about 3 mm should be left behind the knot.

Like all foreign bodies, polyamide threads can give rise to calculus formation when in prolonged contact with the biliary or urinary tract.

When polyamide threads (in particular NYLON) are used in ophthalmic surgery, the following warnings should be borne in mind: patients undergoing eye operations require prolonged observation; corneal polyamide threads should be removed after 3 to 12 months; laser irradiation accelerates the tissue breakdown of polyamide threads.

Monofilament threads (SERALON®, NYLON) are preferable to multifilament threads for use as sutures in infected tissue. A multifilament thread structure can prolong infection.

Surgical instruments should be handled with care in order to avoid damage to threads.

Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care.

### ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of polyamide threads: local irritation, inflammatory reactions to the foreign body; formation of foreign-body giant cells, suture fistulas or granulomas.

### STERILITY

The materials are delivered sterile (sterilization method: ethylene oxide; for SUPRAMID in multipacks: gamma irradiation). Do not resterilize, critical changes in suitability for use possible (more details available on request): Do not use damaged packs!

After opening cassette packs aseptic removal is guaranteed for a period of 6 months. Please make a note of the use by date using the attached pack label.

Packs that have been opened for use but then not used — or in the case of the cassette packs not handled correctly (see INSTRUCTIONS FOR REMOVAL FROM THE PACK) — are to be discarded!

### STORAGE CONDITIONS

No special storage conditions are required. Do not use after the expiry date or, with cassette packs, after the use by date!

### DISPOSAL INFORMATION

Dispose of unused products or waste material with particular care due to the potential risk of contamination

and/or injury. In all other cases observe the specific national regulations.

### INSTRUCTIONS FOR REMOVAL FROM THE PACK

**Minibox:** For easy removal, the thread should be withdrawn from the Minibox not obliquely, but vertically upwards.

**Minibox - microsutures in protective sleeve:** Remove the sleeve from the minibox attachment; grasp the needle with the needle holder and withdraw the thread from the protective sleeve; if a **short thread** is required, it should be cut through together with the protective sleeve; in the case of **double-armed threads**, two individual threads can, if required, be withdrawn after the loop of thread at the end of the protective sleeve has been cut through.


**Cassette pack (also called suture dispenser or outer pack) :** Remove the tamperproof closure (dispenser cap closure) and use the enclosed label to note the expiry date on the outer pack; observe the following instructions to guarantee aseptic removal for a period of 6 months after opening. Flip up the dispenser cap at the clip and fix it to the back part of the cap. If the outer pack has already been opened, i.e. the dispenser cap closure has been removed or opened, **disinfect** the inner surface of the dispenser **before use** with a commercially available alcohol-based product. Pull out the end of the suture with a sterile instrument. Discard the first 10 cm of the suture. Pull out the desired length of suture, cut it off flush, and close the cap. **After finishing the procedure, but no more than 4 hours later, or if required again disinfect** the inner surface of the dispenser and the cap with a commercially available alcohol-based product and close the cap.

### PRESENTATION


Polyamide suture materials are supplied in various thread structures, gauge sizes, lengths, and colours and with various types of needle or without needles. The combinations are supplied in individual packs and multiple packs in a sterile package. Economical cassette packs containing SERALON® und SUPRAMID are available for bulk purchasers. More precise details are given in the catalogue.

## EXPLANATION OF SYMBOLS ON THE PACKAGE

 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

 Surgical suture, nonabsorbable, monofilament, undyed

 Surgical suture, nonabsorbable, monofilament, dyed

 Surgical suture, nonabsorbable, twisted, coated, undyed

 Surgical suture, nonabsorbable, twisted, coated, dyed

 Polyamide

 Antireflective needle

 Detachable needle

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG