**SERAPID®**

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**POLYGLYCOLIC ACID SUTURE MATERIAL****• SERAPID®****DESCRIPTION**

SERAPID® is a braided, synthetic, absorbable, sterile surgical suture material. It is made from the homopolymer of glycolic acid and has the molecular formula  $[C_2O_2H_2]_n$ . It is supplied undyed. SERAPID® is more rapidly absorbed than Serafit® due to the use of polymers with a lower molecular weight.

The suture material complies with the underlying requirements of the Medical Devices Directive 93/42/EEC and the European Pharmacopoeia (EP) monograph "Sutures, sterile synthetic absorbable braided (Fila resorbilia synthetica torta sterilia)".

**• SERAPID®**

polyglycolic acid, braided, coated with calcium stearate and polycaprolactone or polyglycolic acid/poly lactide

**INDICATIONS**

SERAPID® is intended for use in adaptation of soft tissues where short-term absorbable thread material is required.

**USE**

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the doctor, the surgical technique, and the type and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security. In order to prevent needle damage during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint.

**ACTION**

The function of a medical suture in wound care is to attach tissue to tissue.

SERAPID® is degraded by hydrolysis in tissues to form glycolic acid, which is then metabolized in the body. This absorption becomes manifest firstly as a reduction in tensile strength and later as a loss of mass of the thread. The breaking load falls to 50% of its initial value after about 5-7 days. Degradation is complete after 42 days.

**CONTRAINDICATIONS**

SERAPID® must not be used where more prolonged or long-term stability of the suture material is required. SERAPID® is not indicated for ligatures or for use in cardiovascular or surgery or neurosurgery.

**WARNINGS / PRECAUTIONS / INTERACTIONS**

The requirements of the suture material and the techniques required vary with the indication. Therefore, the user should be familiar with the surgical techniques to be used. In the choice of suture material, account should be taken of the in-vivo characteristics of the suture material concerned. In patients with delayed wound healing SERAPID® should be used only after careful consideration. Like all foreign bodies, SERAPID® can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract). Additional use of non-absorbable suture materials or suitable immobilization techniques may be required with sutures that are under tension or require greater support. Use of suitable ancillary measures is required with sutures in infected tissue, as the multifilament structure of the thread can prolong infection. Because of the increased risk of infection associated with the multifilament structure of the thread, skin sutures should if possible be intracutaneous and the suture material, in particular the knots, should be placed as deep as possible. Surgical instruments should be handled with care in order to avoid damage to the thread. Bending of needles or grasping them outside of the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled and disposed of with particular care.

## ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of SERAPID®: transient local irritation, transient inflammatory reactions to the foreign body; rarely, formation of suture fistulas or granulomas, formation of calculus after prolonged contact with salt-containing solutions.

## STERILITY

SERAPID® is supplied sterile (sterilization method: gamma irradiation). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packages. Packages that have been opened for use but then not used are to be discarded.

## STORAGE CONDITIONS

No special storage conditions are required. Do not use after the expiry date.

## DISPOSAL INFORMATION

Dispose of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

## PRESENTATION

Suture materials of SERAPID® are supplied in various gauge sizes and lengths and with various types of needle or without needles. The threads or combinations are supplied in individual packs and multiple packs in a sterile package. More precise details are given in the catalogue.

## EXPLANATION OF SYMBOLS ON THE PACKAGE

 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

Surgical suture, absorbable, braided, coated, undyed

Polyglycolic acid

Antireflective needle

Detachable needle

Needle identification symbols and abbreviations are explained in the catalogue.

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SERAG-WIESSNER GmbH & Co. KG