CE 1014

ce100011/N01935v2

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SILICONE-LOOPS

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).
- The product is intended for single use only.
- Do not resterilize or re-use the product as critical changes in performance and safety are possible.
- Do not use after expiration date.
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or accidentally opened prior to use.
- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

DESCRIPTION AND INTENDED PURPOSE

Silicone Loops are flat, elastic bands. They are made of silicone rubber. Silicone Loops are available in lengths of 45 and 70 cm and in band widths from 1.5 to 5 mm as well as in various colours (White, red, blue or yellow).

Silicone Loops are intended for labelling, encircling or temporary closure of selected tissue structures such as blood or nerve vessels, ureter or vas deferens.

CONTRAINDICATIONS

The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

Silicone Loops are removed before the end of the surgical procedure. They are not intended for implantation. The device is not intended specifically for use in direct contact with the heart, the central circulatory system or the central nervous system.

ADVERSE REACTIONS

None known.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG