de GEBRAUCHSANWEISUNG**NAHTMATERIAL aus POLYESTER in der Großpackung****en INSTRUCTIONS FOR USE****SUTURE MATERIALS of POLYESTER in cassette packs****fr MODE D'EMPLOI****MATERIAU DE SUTURE en POLYESTER en cassette****it ISTRUZIONI PER L'USO****MATERIALE DA SUTURA in POLIESTERE in confezione grande****es INSTRUCCIONES DE USO****MATERIAL DE SUTURA de POLIESTER en el envase clínico****nl GEBRUIKSAANWIJZING****CHIRURGISCH HECHTMATERIAAL uit POLYESTER in de grootverpakking****cs NÁVOD K POUŽITÍ****CHIRURGICKÝ ŠICÍ MATERIÁL z POLYESTERU ve velkém balení****ru ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ****ШОВНЫЕ МАТЕРИАЛЫ ИЗ ПОЛИЭФИРНЫХ ВОЛОКОН в большой упаковке****mk УПАТСТВО ЗА УПОТРЕБА****МАТЕРИЈАЛИ ЗА ШИЕЊЕ ОД ПОЛИЕСТЕР во големо пакување****hu HASZNÁLATI UTASÍTÁS****POLIÉSZTERBŐL KÉSZÜLT VARRÓANYAGOK nagy kiserelésben****sr UPUTSTVO ZA UPOTREBU****MATERIJAL ZA ŠIVENJE OD POLIESTERA u velikom pakovanju****hr UPUTE ZA UPORABU****ŠIVAĆI MATERIJAL od poliestera u skupnoj ambalaži**

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POLYESTER SUTURE MATERIALS in cassette packs

- **SULENE®**
- **TERYLENE**

WARNINGS AND PRECAUTIONS

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: gamma irradiation). 
- The product is intended for single use only. 
- Do not resterilize or re-use the product as critical changes in performance and safety are possible. 
- Do not use after expiration date. 
- For cassette packs, the shelf life after opening should be considered.
- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use. 
- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

DESCRIPTION AND INTENDED PURPOSE

SULENE®, and TERYLENE are synthetic, braided, non-absorbable surgical suture materials of polyethylene terephthalate (PET). They are supplied in an isopropanol-containing preserving solution.

The suture materials are coated with silicone. They are supplied dyed green with D&C No. 6, colour index no. 61565, or undyed. They are available in various gauge sizes and lengths. The suture materials are only supplied without needles. They comply with the current edition of the European Pharmacopoeia monograph (EP) for "Sutures, sterile non-absorbable" (Fila non resorbilia sterilia)". In deviation from this standard, the maximum thread length in case of the cassette packs is not 4 but 100 m.

SULENE®, and TERYLENE are suitable for adaptation of soft tissues and for ligatures where non-absorbable suture material is indicated, in particular where long-term stability of the suture material is required. Polyester threads can also be used as holding sutures and for marking purposes.

INSTRUCTIONS FOR USE

In the choice and use of suture material, due account should be taken of the condition of the patient, the experience of the physician, the surgical technique, and the type and size of the wound. Standard knot-tying techniques should be used to ensure knot-holding security.

Unneeded suture material may be used for ligation or may be connected with surgical eye needles. For a list of suitable needles and application instructions, please refer to

the Instructions for Use of Serag-Wiessner sterile eye needles. The suitability of other manufacturers' products has not been tested. When using other brands, their specific instructions must be observed.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to foreign material. Ligatures are a special type of suture used to close hollow organs. In the body, PET suture material is slowly and progressively encapsulated by connective tissue. It is non-absorbable and does not lose tensile strength in tissues even after prolonged periods.

CONTRAINDICATIONS

SULENE[®], and TERYLENE in the cassette pack are not intended for use in the central vascular or nervous systems. The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques needed vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. Like all foreign bodies, polyester threads can give rise to calculi when in prolonged contact with the biliary or urinary tract. In the case of sutures in infected tissue the use of suitable ancillary measures is called for, as the multifilament thread structure can prolong infection. Surgical instruments should be handled with care in order to avoid damage to threads.

ADVERSE REACTIONS

The following adverse reactions can occur in association with correct use of the device: local irritation, inflammatory reactions to the foreign body; formation of foreign-body giant cells, suture fistulas or granulomas.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

INSTRUCTIONS FOR REMOVAL FROM THE PACK

Remove the tamperproof closure (dispenser cap closure) and use the enclosed label to note the expiry date on the outer pack; observe the following instructions to guarantee aseptic removal for a period of 6 months after opening. Flip up the dispenser cap at the clip and fix it to the back part of the cap. If the outer pack has already been opened, i.e. the dispenser cap closure has been removed or opened, **disinfect** the inner surface of the dispenser **before use** with a commercially available alcohol-based product. Pull out the end of the suture with a sterile instrument. Discard the first 10 cm of the suture. Pull out the desired length of suture, cut it off flush, and close the cap. **After finishing the procedure, but no more than 4 hours later, or if required again disinfect** the inner surface of the dispenser and the cap with a commercially available alcohol-based product and close the cap.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regu-

lations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG