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**KNOCHENWACHS**

**en INSTRUCTIONS FOR USE**

**BONE WAX**

**fr MODE D' EMPLOI**

**SUBSTITUT OSSEUX EN PATE**

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**CERA PER OSSA**

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**WOSK KOSTNY**

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**KOSTNÍ VOSK**

**tr KULLANIM TALİMATI**

**KEMİK BALMUMU**

**lt NAUDOJIMO INSTRUKCIJA**

**KAULŲ VAŠKAS**

**lv LIETOŠANAS INSTRUKCIJA**

**KAULU VASKS**

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**CEARĂ OS**

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## BONE WAX

### DESCRIPTION

Bone wax is a sterile, non-resorbable, soft pliable mass with the following composition:

white beeswax BP 89 % w/w

isopropyl palmitate BP 11 % w/w

Bone wax is supplied sterile in rectangular slabs.

### INDICATIONS

Bone wax is designed to control local bleeding from raw bone surfaces.

### USE

Bone wax is applied sparingly to the affected bony structures. Kneading the wax under sterile conditions will warm it until the desired consistency is obtained.

### ACTION

Bone wax prevents bleeding from raw bone surfaces by forming a mechanical barrier.

### CONTRAINDICATIONS

Bone wax should not be used in the presence or likelihood of infection in the area to be treated or where rapid bone regeneration is desired.

### WARNINGS / PRECAUTIONS / INTERACTIONS

Bone wax should be prepared and used under strictly aseptic conditions. It should be applied sparingly and excess material removed. Its barrier effect may impair the process of bone healing and obstruct the treatment of any concurrent infection.

### ADVERSE REACTIONS

Bone wax may cause a mild inflammatory reaction in the surrounding tissue and inhibit osteogenesis in the affected area.

### STERILITY

Bone wax is supplied sterile (sterilization method:gamma irradiation). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Packages that have been opened for use but then not used are to be discarded.

### STORAGE CONDITIONS

The storage conditions apply as indicated on the label. It should not be used after the expiry date.

### DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

### PRESENTATION

Bone wax is supplied sterile in packs of 24 x 2.5 gram slabs.

## **EXPLANATION OF SYMBOLS**

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG