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 Monofilament, resorbierbar, gefärbt.

**Synthetic surgical textile implant
made of POLYDIOXANONE**
• SERASYNTH® MESH

DESCRIPTION

SERASYNTH® MESH is an absorbable sterile surgical implant made from a synthetic monofilament ground thread. This thread is composed of the polymer poly-p-dioxanone with the molecular formula $[C_4H_6O_3]_n$ and is available dyed with D+C violet No. 2, Colour Index No. 60725. The thread material used conforms to SERASYNTH® suture material.

INDICATIONS

SERASYNTH® MESH is intended for use as temporary surgical support, bridging and/or strengthening bodily structures, when a long-term (up to six weeks) absorbable textile implant is indicated.

USE

Different dimensions are available for the various indications. Furthermore, the mesh can be cut-to-size to individually adapt to different anatomical conditions.

ACTION

SERASYNTH® MESH strengthens and bridges temporarily damaged body structures. It is degraded in the tissue by means of hydrolysis to 2-hydroxyethoxy acetic acid, which is subsequently metabolised by the organism. This absorption expresses itself initially as a loss of tensile strength of the mesh, which is later accompanied by a loss of mass. After two weeks, tensile strength is reduced to approximately 70-80% and after four weeks to about 50-70% of the original amount. Absorption is complete after about 30 weeks.

CONTRAINDICATIONS

SERASYNTH® MESH should not be used where a lasting stability of the implant is required.

WARNINGS / PRECAUTIONS / INTERACTIONS

The user should be familiar with the various surgical techniques and procedures. When selecting the implant, its in-vivo properties should be taken into account. The use of SERASYNTH® MESH should be carefully considered in patients with delayed wound healing.

Implants should always be applied under consideration of the local anatomy, thereby taking particular care to avoid damaging larger blood vessels as well as nerves and organs.

Implants that are under tension or those needing additional support may require the supplemental use of non-absorbable materials or suitable immobilisation procedures.

The surgical instruments should be handled carefully to avoid damage to the textile implant or to any suture materials used for fixation. Textile implants can, for example, be mechanically damaged by staples, clips or clamps. The textile implants may not be fixated in the outermost mesh on the edge of the implant; a distance of at least 6 mm or 2 meshes to the mesh edge must be maintained.

ADVERSE REACTIONS

The following adverse reactions may occur even with the correct use of SERASYNTH® MESH: Transient tissue reactions in the form of localised irritation, inflammatory foreign-body reactions and seroma formation.

STERILITY

SERASYNTH® MESH is supplied in a sterile pack (sterilisation method: ethylene oxide). Do not resterilize, critical changes in suitability for use possible (more details available on request)! Do not use damaged packs! Discard unused opened packs!

STORAGE CONDITIONS

No special storage conditions are required. Do not use after the expiry date.

DISPOSAL INFORMATION

Dispose of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

PRESENTATION

SERASYNTH® MESH is available in various sizes and dimensions. If required, it can be provided with incisions for easy adaptation to anatomical structures. Dimensions < 15 x 15 cm are packed flat, while the other sizes are packed folded. Please consult our catalogue for more information.

EXPLANATION OF SYMBOLS ON THE PACKAGE

 CE 1014 CE symbol and identification number of the notified body. The product satisfies the essential requirements of Directive 93/42/EEC on medical devices

 PDO Polydioxanone

 Textile mesh implant with different structures and dimensions.

The symbol within the circle shows the basic structure of the mesh implant.

In detail this can be:

 Monofilament, absorbable, dyed

 Do not use if package is damaged.

 Do not re-use

 STERILE EO Sterilized with ethylene oxide

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SERAG-WIESSNER GmbH & Co. KG