SERASIS[®] - / SERAPRO[®] - INSTRUMENT

de

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Reusable instruments made of surgical steel SERASIS® and SERAPRO® Instruments

DESCRIPTION

SERASIS® and SERAPRO® instruments are reusable instruments made of surgical steel.

INDICATIONS

The various reusable instruments allow insertion of SERASIS®, SERATOM®, SERATEX® or SERATIM® implants, which are used in incontinence and prolapse surgery.

Depending on the procedure and indication, there are various instrument modifications to choose from (see product options).

USE

For instructions on assembly and use, please refer to more detailed information, which describes how to handle the implantable products.

WARNINGS / ADVERSE REACTIONS

The instruments should only be used by people who have had training in instrument handling. Under no circumstances may the instruments made of hardened steel be bent into shape or reworked, as there is a risk they may otherwise break. Dispose of bent or reworked instruments.

The passage of instruments can damage blood vessels, nerves, the bladder or the bowel in the form of punctures or tears that require surgical repair.

CLEANING, DISINFECTION, STERILISATION and PREPARATION FOR REUSE

The instruments are supplied unsterile and must be cleaned, disinfected and sterilised before they are used for the first time and before every subsequent use.

When preparing medical devices, the following steps should be taken if required: appropriate preparation, cleaning / disinfection / rinsing / drying, checking for clean / intact surfaces, maintenance / repair, checking safety functions, labelling, packaging and sterilisation.

Cleaning and disinfection can be automatic or manual and include or exclude ultrasound treatment. The chemicals to be used may be acidic, neutral or alkaline with or without tensides. The temperature must not exceed 60°C when using chemicals and 93°C when using DI water (thermal disinfection). It is assumed that commercially available cleaning and disinfection products that have been approved for the application are used in the recommended concentrations, for the recommended periods of time and at the recommended temperatures. The instruments must be free of residues (particularly protein). A final rinse with DI water followed by drying at no more than 100°C is recommended.

The following recommendations apply to sterilisation:

Steam sterilisation: as described in the European Pharmacopoeia (121°C, 20 minutes) or other combinations of holding time and temperature (e.g. 134°C, 18 minutes) that also guarantee a sterility assurance level (SAL) of at least 10^{-6} , whereby the specified conditions are maximum values.

Other methods: formaldehyde sterilisation (max. 65°C), EO sterilisation (max. 65°C), gas plasma sterilisation (e.g. STERRAD®) or liquid sterilisation (e.g. STERIS®), which have been validated with an SAL of at least 10^{-6} by the user.

Preparation for reuse:

After use, the instruments may be kept either wet (immediately after the operation, dismantle the instruments, if required, and lay them in a non-fixing, active cleaning disinfection solution) or dry (without disinfectant or other additional liquids), prior to sending them for cleaning, disinfection and sterilisation (see above).

Before use, the medical devices should be examined for any loss of function or visual changes that prohibit their reuse.

Irrespective of the sterilisation conditions described above, reuse of medical devices should be ruled out and the products discarded if there is any suspicion of prion contamination (see Robert Koch Institute recommendation).

The use of cleaners containing phosphoric or hydrofluoric acid can cause laser marks on the products to fade or disappear.

The information given here cannot replace detailed process descriptions, as it is impossible to go into detail about the many treatment methods used across the world. The above information is supplied without liability and does not absolve users of their responsibility.

STORAGE CONDITIONS

The reusable instruments should be stored dry, clean, and in the unclamped position.

PRESENTATION

SERASIS® and SERAPRO® instruments are supplied as unsterile reusable products and available in several modifications for different indications. For further details, please see our catalogue.

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SERAG-WIESSNER GmbH & Co. KG