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POLYGLYCOLIC ACID-CAPROLACTONE SUTURE MATERIAL**● SERAFAST®****WARNINGS AND PRECAUTIONS**

- Read the instructions for use in their entirety before using the product.
- The product should only be used by qualified and trained personnel. It is the physician's responsibility to be familiar with both appropriate surgical techniques and anatomy prior to use.
- The product is delivered sterile (sterilization method: ethylene oxide).

- The product is intended for single use only.

- Do not resterilize or re-use the product as critical changes in performance and safety are possible.

- Do not use after expiration date.

- Do not use if the product shows signs of damage (e.g. cuts or kinks) or if the outer packaging or sterile packaging has been damaged or unintentionally opened prior to use.

- Use the product only for the intended indication.
- Please contact the manufacturer for further information.

DESCRIPTION AND INTENDED PURPOSE

SERAFAST® is a synthetic absorbable monofilament surgical suture material. It is composed of a segmented copolymer with the basic compounds glycolide and ϵ -caprolactone. It is available dyed with D+C violett No. 2, Colour Index No. 60725, or undyed. SERAFAST® is available in various gauge sizes and lengths. It is supplied with various types of needles made of surgical steel. SERAFAST® may be supplied in combination with other components, such as a suture guide tube.

The suture material complies with the underlying requirements of the European Pharmacopoeia (EP) monograph "Sutures, sterile synthetic absorbable monofilament (Fila resorbilia synthetica monofilamenta sterilia)".

SERAFAST® is used to approximate soft tissues or for tying ligatures in cases where rapidly absorbable (up to 2 weeks) suture is indicated.

INSTRUCTIONS FOR USE

In the choice and use of the suture material, due account should be taken of the condition of the patient, the experience of the physician, the surgical technique, and the type and size of the wound. Knot security requires the established standard knotting

techniques for synthetic monofilament sutures. In order to prevent needle damage during suturing it is recommended that the needle be grasped at an area between the end of its first third (as seen from its end) and its midpoint.

ACTION

The function of a medical suture in wound care is to attach tissue to tissue or tissue to foreign material. Ligatures are a special type of suture intended to close hollow organs. SERAFAST® is broken down in the tissue by hydrolysis to adipic acid which is then metabolised in the body. This absorption expresses itself as a loss of tensile strength which is later accompanied by loss of mass. The breaking strength is reduced to about 60% of its initial value within 7 days (in-vitro). Absorption is completed after about 90 days (in-vivo).

CONTRAINDICATIONS

SERAFAST® must not be used if longer-term or permanent stability of the suture material is required, e.g. securing synthetic implants (heart valves etc.). SERAFAST® is not indicated for use in the central circulatory or central nervous system.

The device must not be used in patients with allergies to any of the product components.

RESTRICTIONS AND LIMITATIONS OF USE / NOTES / INTERACTIONS

The requirements of the suture material and the techniques required vary with the indication. In the choice of suture material, account should be taken of the characteristics of the suture material concerned. The use of SERAFAST® should be carefully considered in patients with delayed wound healing. Like all foreign bodies, the device can give rise to calculus formation when in prolonged contact with salt-containing solutions (e.g. in the biliary or urinary tract). Sutures under tension or those needing additional support may require the supplemental use of non-absorbable sutures or suitable immobilisation procedures.

In intracutaneous sutures the suture material – in particular the knots – should be placed as deep as possible. Surgical instruments should be handled with care in order to avoid damage to the thread. Bending of needles or grasping them outside the recommended area can result in needle breakage. Because of the risk of infection from needle pricks, needles should be handled with particular care.

Considerations should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur. The use of absorbable sutures may be inappropriate in elderly, malnourished or debilitated patients. Skin sutures which remain in place longer than 7 days may cause local irritation and should be snipped off or removed.

ADVERSE REACTIONS

The following adverse reactions may occur in association with correct use of the device:

local irritation, inflammatory reactions to the foreign body.

Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

STORAGE CONDITIONS

The storage conditions apply as indicated on the label.

DISPOSAL INFORMATION

Disposal of unused products or waste material with particular care due to the potential risk of contamination and/or injury. In all other cases observe the specific national regulations.

MEANING OF THE SYMBOLS

Symbols are explained at the end of the instructions for use.

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SERAG-WIESSNER GmbH & Co. KG