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---

## **Description**

Serag BSS is a physiological salt solution, adapted to the aqueous humour of the eye.

## **Composition**

1000 ml solution contain:

Sodium chloride	6.400 g
Potassium chloride	0.750 g
Calcium chloride 2 H <sub>2</sub> O	0.480 g
Magnesium chloride 6 H <sub>2</sub> O	0.300 g
Sodium acetate 3 H <sub>2</sub> O	3.900 g
Sodium citrate 2 H <sub>2</sub> O	1.700 g

Hydrochloric acid (for pH adjustment)

Water for injections

## **Therapeutic indications**

- Irrigation of the anterior chamber of the eye during cataract surgery and other intraocular procedures
- Removal of alpha-chymotrypsin after zonulysis
- Intraoperative irrigation of the anterior chamber of the eye in the event of anterior chamber haemorrhage
- Irrigation of the cornea during strabotomy and other extraocular procedures
- Irrigation of the conjunctiva following fluorescein instillation
- Irrigation of the anterior area of the eye and irrigation of the lacrimal passages following chemical burns and burns
- Moistening of the cornea and conjunctiva during laser treatment

---

### Precautions for use and warnings

Do not use for injection or infusion.

In patients with diabetes, the aqueous humour is hyperosmolar and intraoperative lens changes can occur. As Serag BSS is iso-osmotic with the normal aqueous humour, it should be used with caution in diabetic patients undergoing vitrectomy. In this case, it is recommended to align the osmolarity of the aqueous humour.

Do not resterilize. Critical change of usability is possible.

### Interactions with other medicines

Serag BSS is adjusted to the osmotic conditions of the normal aqueous fluid. Do not add any other substances.

### Dosage instructions, type and duration of administration

The dosage depends upon the specific circumstances of the operative, therapeutic or diagnostic procedure and must be determined by the physician on an individual basis.

For irrigation of the eye during surgical procedures and diagnostic and therapeutic measures requiring irrigation or moistening. The duration of use is limited to the corresponding procedure.

### Overdose

The risk of poisoning (intoxication) is excluded with regard to administration of Serag BSS, as the product is a physiological irrigation solution adjusted to the aqueous humour.

### Side effects

When the corneal endothelium is abnormal, any intraocular procedure (e.g. intraocular irrigation) can lead to corneal damage (keratopathy), characterised by corneal oedema or clouding. Post-operative inflammatory reactions can occur. To date, their relation to the use of Serag BSS has not been established.

### Shelf life and storage conditions

Do not use after the expiry date!

Only use clear solutions!

Use immediately after opening the container! Discard any unused product.

### **Disposal informations**

Dispose of unused products or waste material in accordance with the specific national regulations.

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SERAG-WIESSNER GmbH & Co. KG