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LAVANID® - Wundgel

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LAVANID® - Wound Gel

DESCRIPTION AND PRODUCT PROFILE

The surface of chronic wounds is generally made up of devitalized and necrotic tissue, and is also colonised with bacteria. For wound treatment to be successful the wound surface needs to be removed, e.g. through surgical debridement. Wound irrigation and moistening are a necessary part of this process. Chronic wounds exhibit dysfunctional wound healing processes. This is why treatment with preparations that are very well tolerated by the tissues is indicated.

LAVANID® Wound Gel is a sterile, preserved Ringer's-solution-based wound gel which stands out due to its very good tissue tolerance. Pol-hexanide 0.04% is used as a preservative. Polyhexanide is evaluated as the agent of choice for chronic wounds and burns. In addition to cleaning and moistening wounds, LAVANID® Wound Gel is also appropriate for preservative dressings, e.g. for wounds colonised with multiresistant bacteria. Preservative dressings reduce the risk of carrying the bacteria from the wound to the patient's environment and they protect the wound from bacteria getting in from the outside.

COMPOSITION
Ringer's solution (isotonic electrolyte solution made up of sodium chloride, potassium chloride, calcium chloride 2 H₂O, water for injection), Glycerin, macrogol, hydroxyethylcellulose, polyhexanide 0.04%.

The ingredient polyhexanide can cause allergic reactions. There are individual case reports of acute immune system reactions (anaphylactic reactions). Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

TISSUE TOLERANCE
Experimental and clinical studies have shown that the combination of Ringer's solution with 0.04% polyhexanide is very well tolerated by tissue. Clinical trials with LAVANID® Wound Gel have shown that it does not negatively affect the formation of healthy new tissue. There are no local signs of irritant conditions. Use of this product is pain-free for patients.

Because it is so well tolerated by tissues LAVANID® Wound Gel is suitable for repeated and long-term wound treatment.

LAVANID® Wound Gel is sterile as long as the packaging is unopened and undamaged. Although it is no longer sterile after the container has been opened, because it contains a preservative it can be used up to one month after opening under certain conditions. These conditions are that the contents be protected from contamination. The container must not come in contact with the wound. The opening should be capped immediately after each use. Gel residue around the opening should be wiped off with a sterile compress.

INSTRUCTIONS FOR USE

- Prepare and clean the wound as usual.
- Apply LAVANID® Wound Gel to the wound or a wound dressing (e.g. pads or compresses). Dress wound as usual. LAVANID® Wound Gel should remain on the wound until the next dressing change. For optimal wound moistening we recommend increasing the amount of LAVANID® Wound Gel and lengthening the intervals between dressing changes. If the dressing sticks to the wound while being changed, moisten the dressing with a rinsing solution (e.g. LAVANID® Wound Rinsing Solution) before removing it.
- If the LAVANID® Wound Gel is going to be used more than once, please read the shelf-life information.

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- LAVANID® Wound Gel must not be used
 - in cases of known hypersensitivity to any of the ingredients,
 - in the area of cartilage and joints.
- There is insufficient data on its use during pregnancy and breast-feeding. Therefore, LAVANID® Wound Gel should not be used unless urgently indicated.

WARNINGS AND PRECAUTIONS FOR USE
For external use only. Do not mix with additives. Sterile if package is undamaged. Do not resterilize! This may critically alter its usability.

Do not use after the expiry date.

INTERACTIONS WITH OTHER DRUGS
When used in combination with other wound treatment products, the instructions for use of these products should be checked for possible interactions with the ingredients of LAVANID®. Due to incompatibility problems with the ingredient polyhexanide, LAVANID® Wound Gel must not be used in combination with anionic substances (e.g. anionic tensides, anionic gel formers and products containing PVP-iodine). In instances where there is prior contact, these substances must be rinsed liberally out of the wound.

SIDE EFFECTS
The ingredient polyhexanide can cause allergic reactions. There are individual case reports of acute immune system reactions (anaphylactic reactions). Users and/or patients are encouraged to report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority.

INFORMATION ABOUT SHELF-LIFE
The expiry date can be found on the container. Do not use after the expiry date!

INDICATIONS
LAVANID® Wound Gel is used to moisten and clean wounds and to moisten bandages and wound dressings with preservative.

INFORMATION ABOUT SHELF-LIFE
The expiry date can be found on the container. Do not use after the expiry date!

Shelf life after opening: 1 month

LAVANID® Wound Gel is sterile as long as the packaging is unopened and undamaged. Although it is no longer sterile after the container has been opened, because it contains a preservative it can be used up to one month after opening under certain conditions. These conditions are that the contents be protected from contamination. The container must not come in contact with the wound. The opening should be capped immediately after each use. Gel residue around the opening should be wiped off with a sterile compress.

DISPOSAL INFORMATION
Dispose of unused products or waste material in accordance with the specific national regulations.

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