

Safe and gentle spot landing with the

LAVANID® dressing

Cutting edge
wound management



Sterile, moist
dressing with
polihexanide



150
YEARS
1866 - 2016



SERAG
WIESSNER

Building blocks for

systematic wound management



Optimal wound care:

Clean the wound with
LAVANID® products



Primary wound cover
with a LAVANID® dressing



Secondary dressing



Fixation

Optimal wound care

A thorough cleaning of the wound is a basic requirement for optimal wound care and fundamental to rapid healing. LAVANID® dressings now perfectly complement the LAVANID® 1/2 wound irrigation solutions and the LAVANID® / LAVANID® V+ wound gels.

The new product puts the finishing touch to an optimally coordinated system for primary wound care. We can therefore provide a full wound care system from cleaning to the primary covering of acute and chronic wounds.


Description and product profile

The LAVANID® dressing is a ready-to-use sterile, moist dressing with polihexanide.

After thoroughly cleaning the wound with LAVANID® wound irrigation solutions and applying LAVANID® wound gel if required, the LAVANID® dressing is applied as primary covering.

The primary dressing consists of a biocelulose matrix with an antimicrobial complex. This complex in the dressing provides efficient infection prophylaxis against micro-organisms, protecting against bacterial and fungal pathogens.

LAVANID® dressings dissolve biofilms to create an optimal milieu for the formation of granulation tissue in acute and chronic wounds. The wound contact meshes on both sides of the dressing prevent it from adhering to the wound, which makes dressing changes easier and more comfortable.





Ready-to-use moist dressing for direct application on the wound



Can be easily cut, draped or used as a tampon for packing the wound



Can also be additionally moistened with LAVANID® wound gel (with dry wounds)

Indications for use

LAVANID® dressings are used for acute and chronic wounds with a risk of infection.

For example:

- venous leg ulcers
- diabetic foot syndrome
- pressure sores (up to grade 3)
- postoperative wounds
- burns (up to second degree)
- superficial wounds, such as split-skin graft donor sites

Safe and versatile application

- Can be used on both sides
- Very easy to drape or cut
- Good vertical absorption of exudate
- Rapid and efficient wound debridement
- Can be used in all phases of wound healing

Benefits: LAVANID® dressing

Efficient infection prophylaxis

Dissolves the biofilm and cleans the wound

Good vertical absorption of exudate

Easy and safe handling

Can be combined with conventional secondary dressings!

Wound phase	Exudation			
	dry	normal	moderate	severe
Inflammation phase	LAVANID® dressing	LAVANID® dressing	LAVANID® dressing	LAVANID® dressing
	+ LAVANID® wound gel + film dressing or + semi-permeable wound dressing	+ LAVANID® wound gel if needed + film dressing or + foam dressing + fixation	+ absorptive compress + fixation or + superabsorbent dressing + fixation	+ superabsorbent dressing + fixation
Granulation/epithelisation phase	LAVANID® dressing	LAVANID® dressing	LAVANID® dressing	LAVANID® dressing
	+ film dressing or + semi-permeable wound dressing	+ sterile gauze + film dressing or + semi-permeable wound dressing	+ absorptive compress + fixation or + superabsorbent dressing + fixation	+ superabsorbent dressing + fixation

Ordering information:

Art. no.	PIP code	Size	Sales Unit
016310	12388506	10 x 10 cm	10 pieces



HAWEST Research AG,
CH-6312 Steinhausen / Switzerland
www.hawest.com

Distributed by:



**SERAG
WIESSNER**

SERAG-WIESSNER GmbH & Co. KG
Zum Kugelfang 8 - 12
95119 Naila/Germany



+ 49 9282 937-0



+ 49 9282 937-9369

Export Department:



+ 49 9282 937-230



+ 49 9282 937-9785



info@serag-wiessner.de



www.serag-wiessner.de